LISANTE M-58LEO(I)

[PRODUCTO]

LISANTE M-58LEO(I)

[TAMAÑO DEL ENVASE]

1 l x 4

[USO PREVISTO]

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1L×4

[INTENDED USE]

M-58LEO(I) Lyse breaks down red blood cell membrane and cooperates with the M-58LEO(II) Lyse to 4-differentiate WBCs.

[PRINCIPLE]

The M-58LEO(I) is formulated to lyse red blood cells and process white blood cells in the blood cell measurement, amplify the differences among WBC populations, conducting DIFF channel WBC analysis by incorporating laser scatter flow cytometry.

[ACTIVE INGREDIENTS]

Surfactant<35g/L

[STORAGE AND EXPIRATION DATE]

The product will be stable for 2 years when stored at 2°C to 30°C(35°F to 86°F),and the relative humidity must not exceed 90%. The working temperature range of the product is consistent with that of its applicable instruments. The open-vial validity is 60 days.

[APPLICABLE INSTRUMENT]

This product applies to BC-5600 and BC-5800 Auto Hematology Analyzer produced by Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

[SAMPLE REQUIREMENT]

Fresh human whole blood anticoagulated. Do not use contaminated samples.

[MATERIALS REQUIRER BUT NOT PROVIDED]

The following materials are required but not provided with the product: Mindray-manufactured measurement instruments and matched reagents.

[INSTRUCTIONS FOR USE]

Unpack the reagent container. Insert the pickup pipe from the analyzer's cap assembly into the container, making sure the container's cap and the connected assembly connector are in the same color. Turn and secure the cap assembly tightly, and replace the reagent in accordance with the instructions in the Operator's Manual of the analyzer. For detailed usage instruction, please refer to the Operator's Manual of the analyzer.

[CUT-OFF VALUE/REFERENCE INTERVALS]

Not applicable.

[RESULT ELABORATION]

Not applicable.

[LIMITATIONS]

Not applicable.

[PRODUCT SPECIFICATIONS]

1. Appearance: Transparent liquid without deposits, suspended grains or flocks.
2. Blank count results: The blank count results

tested on the matched Auto Hematology Analyzer shall meet the requirements of Table 1.

Table 1 Blank Count Requirements

Parameter	Blank Count Requirements
WBC	≤0.3×10 ⁹ /L

[PRECAUTIONS]

1. For professional use in vitro diagnosis only.
2. Read the package insert carefully before using this product. It shall be used before the expiration date and disposed of properly when expired.
3. Do not use the reagent if it is frozen.
4. If the blank count is abnormal after the reagent is transported, place it still for 24 hours at room temperature before using it.
5. If the reagent is polluted or affected by other factors and becomes abnormal, stop using it and replace it with a normal one.
6. Dispose of the waste, residual and contaminated packing based on local regulations.
7. The following factors can affect the sample analysis: expired or ineffective reagent; reagent polluted by dust in the air; improper disposal of the sample; mixed with or used with reagents produced by other company; mixed use of residual from the old container and the newly-opened; used in condition other than specified.
8. Take the necessary precautions for the use of the product. Do not swallow. Avoid contact with skin and mucous membranes. If you accidentally take the reagent into your mouth, or the reagents accidentally spill on your skin or into your eyes, wash them off with plenty of water and go seek medical treatment if necessary.
9. Disposal of waste liquid and materials should be in accordance with local guidelines.
- 10.The Material Safety Data Sheet (SDS) is available upon request
11. All identified risks have been reduced as far as possible by generally acknowledged state of art, and the overall residual risk is acceptable.
12. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

[GRAPHICAL SYMBOLS]

LOT Batch code	Use-by date	Temperature limit
HUMIDITY LIMITATION	Manufacturer	In vitro diagnostic medical device
CATALOGUE NUMBER	Authorized representative in the European Community	Consult instructions for use

European Conformity

Unique Device Identifier

[REFERENCES]

Not applicable.

[COMPANY CONTACT]

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, Shenzhen 518057, P.R.China Website: www.mindray.com

E-mail Address: service@mindray.com Tel: +86 755 81888998

Fax: +86 755 26582680

EC-Representative: Shanghai International Holding Corp. GmbH(Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

Tel: 0049-40-2513175

Fax: 0049-40-255726

© 2021 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved.

[APPROVAL DATE OF THE OPERATOR ' S MANUAL]

2021.08.30

M-58LEO(II) LYSE

[PRODUCT]

M-58LEO(II) LYSE

[PACKAGE SIZE]

500mL×4

[INTENDED USE]

M-58LEO(II) Lyse cooperates with the M-58LEO(I) Lyse to 4-differentiate WBCs.

[PRINCIPLE]

The M-58LEO(II) LYSE is formulated to lyse red blood cells and process white blood cells in the blood cell measurement, amplify the differences among WBC populations, conducting DIFF channel WBC analysis by incorporating laser scatter flow cytometry.

[ACTIVE INGREDIENTS]

Surfactant.....<40g/L
Dye.....0.01-0.1g/L

[STORAGE AND EXPIRATION DATE]

The product will be stable for 2 years when stored at 2°C to 30°C(35°F to 86°F),and the relative humidity must not exceed 90%. The working temperature range of the product is consistent with that of its applicable instruments. The open-vial validity is 60 days.

[APPLICABLE INSTRUMENT]

This product applies to BC-5600 and BC-5800 Auto Hematology Analyzer produced by Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

[SAMPLE REQUIREMENT]

Fresh human whole blood anticoagulated. Do not use contaminated samples.

[MATERIALS REQUIRER BUT NOT PROVIDED]

The following materials are required but not provided with the product: Mindray-manufactured measurement instruments and matched reagents.

[INSTRUCTIONS FOR USE]

Unpack the reagent container. Insert the pickup pipe from the analyzer's cap assembly into the container, making sure the container's cap and the connected assembly connector are in the same color. Turn and secure the cap assembly tightly, and replace the reagent in accordance with the instructions in the Operator's Manual of the analyzer. For detailed usage instruction, please refer to the Operator's Manual of the analyzer.

[CUT-OFF VALUE/REFERENCE INTERVALS]

Not applicable.

[RESULT ELABORATION]

Not applicable.

[LIMITATIONS]

Not applicable.

[PRODUCT SPECIFICATIONS]

1. Appearance: Transparent liquid without deposits, suspended grains or flocks.
2. Blank count results: The blank count results tested on the matched Auto Hematology Analyzer shall meet the requirements of Table 1.

Table 1 Blank Count Requirements

Parameter	Blank Count Requirements
WBC	≤0.3×10 ⁹ /L

[PRECAUTIONS]

1. For professional use in vitro diagnosis only.
2. Read the package insert carefully before using this product. It shall be used before the expiration date and disposed of properly when expired.
3. Do not use the reagent if it is frozen.
4. If the blank count is abnormal after the reagent is transported, place it still for 24 hours at room temperature before using it.
5. If the reagent is polluted or affected by other factors and becomes abnormal, stop using it and replace it with a normal one.
6. Dispose of the waste, residual and contaminated packing based on local regulations.
7. The following factors can affect the sample analysis: expired or ineffective reagent; reagent polluted by dust in the air; improper disposal of the sample; mixed with or used with reagents produced by other company; mixed use of residual from the old container and the newly-opened; used in condition other than specified.
8. Take the necessary precautions for the use of the product. Do not swallow. Avoid contact with skin and mucous membranes. If you accidentally take the reagent into your mouth, or the reagents accidentally spill on your skin or into your eyes, wash them off with plenty of water and go seek medical treatment if necessary.

9. Disposal of waste liquid and materials should be in accordance with local guidelines.

10.The Material Safety Data Sheet (SDS) is available upon request

11. All identified risks have been reduced as far as possible by generally acknowledged state of art, and the overall residual risk is acceptable.

12. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

[GRAPHICAL SYMBOLS]

LOT Batch code	Use-by date	Temperature limit
HUMIDITY LIMITATION	Manufacturer	In vitro diagnostic medical device
CATALOGUE NUMBER	Authorized representative in the European Community	Consult instructions for use
European Conformity	Unique Device Identifier	

[REFERENCES]

Not applicable.

[COMPANY CONTACT]

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, Shenzhen 518057, P.R.China Website: www.mindray.com

E-mail Address: service@mindray.com Tel: +86 755 81888998

Fax: +86 755 26582680

EC-Representative: Shanghai International Holding Corp. GmbH(Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

Tel: 0049-40-2513175

Fax: 0049-40-255726

© 2021 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved.

[APPROVAL DATE OF THE OPERATOR ' S MANUAL]

2021.08.30