LISANTE M-58LH

[PRODUCTO]

LISANTE M-58LH

[TAMAÑO DEL ENVASE]

500 ml × 4

[USO PREVISTO]

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500mLx4

[INTENDED USE]

M-58LH Lyse breaks down red blood cell membrane and converts hemoglobin to a hemoglobin complex to determine the HGB.

[PRINCIPLE]

M-58LH lyse dissolves RBCs and reacts with hemoglobin to form compound for WBC count, WBC DIFF and hemoglobin measurement.

[ACTIVE INGREDIENTS]

Surfactant<50g/L
Sodium Sulfate Anhydrous.....5-10g/L
Glycerol.....2-10g/L

[STORAGE AND EXPIRATION DATE]

The product will be stable for 2 years when stored at 2°C to 30°C (35°F to 86°F), and the relative humidity must not exceed 90%. The working temperature range of the product is consistent with that of its applicable instruments. The open-vial validity is 60 days.

[APPLICABLE INSTRUMENT]

This product applies to BC-5600 and BC-5800 Auto Hematology Analyzer produced by Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

[SAMPLE REQUIREMENT]

Fresh human whole blood anticoagulated. Do not use contaminated samples.

[MATERIALS REQUIRED BUT NOT PROVIDED]

The following materials are required but not provided with the product: Mindray-manufactured measurement instruments and matched reagents.

[INSTRUCTIONS FOR USE]

Unpack the reagent container. Insert the pickup pipe from the analyzer's cap assembly into the container, making sure the container's cap and the connected assembly connector are in the same color. Turn and secure the cap assembly tightly, and replace the reagent in accordance with the instructions in the Operator's Manual of the analyzer.

For detailed usage instruction, please refer to the Operator's Manual of the analyzer.

[CUT-OFF VALUE/REFERENCE INTERVALS]

Not applicable.

[RESULT ELABORATION]

Not applicable.

[LIMITATIONS]

Not applicable.

[PRODUCT SPECIFICATIONS]

- Appearance: Transparent liquid without deposits, suspended grains or flocks.
- Blank count results: The blank count results tested on the matched Auto Hematology Analyzer shall meet the requirements of Table 1.

Table 1 Blank Count Requirements

Parameter	Blank Count Requirements
WBC	≤0.3×10 ⁹ /L
HGB	≤1g/L

[PRECAUTIONS]

- For professional use in vitro diagnosis only.
- Read the package insert carefully before using this product. It shall be used before the expiration date and disposed of properly when expired.
- Do not use the reagent if it is frozen.
- If the blank count is abnormal after the reagent is transported, place it still for 24 hours at room temperature before using it.
- If the reagent is polluted or affected by other factors and becomes abnormal, stop using it and replace it with a normal one.
- Dispose of the waste, residual and contaminated packing based on local regulations.
- The following factors can affect the sample analysis: expired or ineffective reagent; reagent polluted by dust in the air; improper disposal of the sample; mixed with or used with reagents produced by other company; mixed use of residual from the old container and the newly-opened; used in condition other than specified.
- Take the necessary precautions for the use of the product. Do not swallow. Avoid contact with skin and mucous membranes. If you accidentally take the reagent into your mouth, or the reagents accidentally spill on your skin or into your eyes, wash them off with plenty of water and go seek medical treatment if necessary.
- Disposal of waste liquid and materials should be in accordance with local guidelines.
- The Material Safety Data Sheet (SDS) is available upon request
- All identified risks have been reduced as far as possible by generally acknowledged state of art, and the overall residual risk is acceptable.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.
- This product contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

Warning	
H411	Toxic to aquatic life with long lasting effects.
H315	Causes skin irritation.
H319	Causes serious eye irritation.
Prevention	
P273	Avoid release to the environment.
P280	Wear protective gloves, protective clothing, eye protection and face protection.
P264	Wash all exposed external body areas thoroughly after handling.
Response	
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313	If eye irritation persists: Get medical advice/attention.
P391	Collect spillage.
P302+P352	IF ON SKIN: Wash with plenty of water.
P332+P313	If skin irritation occurs: Get medical advice/attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents/container to authorised hazardous or special waste collection point in accordance with any local regulation.

[GRAPHICAL SYMBOLS]

Batch code	Use-by date	Temperature limit
In vitro diagnostic medical device	Manufacturer	Consult instructions for use

Authorized representative in the European Community	Catalogue number	Humidity limitation
European Conformity	Unique device identifier	

[REFERENCES]

Not applicable.

[COMPANY CONTACT]

Manufacturer Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Address Mindray Building Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R. China.
Website www.mindray.com
E-mail address service@mindray.com
Tel +86 755 81888998
Fax +86 755 26582680
EC-Representative Shanghai International Holding Corp. GmbH(Europe)
Address EiffeustraÙe 80, 20537 Hamburg, Germany
Tel 0049-40-2513175
Fax 0049-40-255726

© 2021-2022 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved.

[APPROVAL DATE OF THE OPERATOR'S MANUAL]

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© 2021-2022 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Tüm Hakları Saklıdır.